

Biopsy Needle Guides for Mindray

- Integral design without any accessories
- High grade medical stainless steel
- Artificial polishing
- Supporting accurate and angle adjustable biopsy
- Matching a wide range of needle gauges
- Simple pressing, quick release



About Us

As a China national high-tech enterprise, InnoFine are dedicated to R&D, manufacturing and sales ultrasound interventional products and consumables medical products since its establishment in 2013, providing safe and reliable, cost-effective products and overall solutions for thousands of hospitals in more than 100 countries around the world, products covering Ultrasound Needle Guides, Probe Cover, Gel, Biopsy Needle etc.

InnoFine will continue to adhere to the spirit of "Be innovative and forge ahead, Strive for excellence and perfection, Customer oriented ", and make unremitting efforts to promote the development of global interventional ultrasound field.

Cooperation Ultrasound Brand (Part)

Reference Sites (Part)

 Mount Elizabeth Novena Hospital,Singapore	 Apollo Hospitals,India	 Semmelweis University Medical Center,Hungary	 Kinderwunschzentrum Heidelberg,Germany	 Orchid Fertility,UAE
 Chinese PLA General Hospital	 Peking Union Medical College Hospital	 West China Hospital Sichuan University	 Shengjing Hospital China Medical University	 Tongji Hospital,Huazhong University of Science and Technology
 Shanghai Ninth People's Hospital	 Xiangya Hospital Central South University	 The First Affiliated Hospital Sun Yat-sen University	 Shenzhen People's Hospital	 Shenzhen Hospital Peking University

Classial Product (Part)



Biopsy needle guides for Mindray SC6-1U probe (Convex type probe)




Biopsy needle guides for Mindray SP5-1U probe (Phased array type probe)



Biopsy needle guides for Mindray V11-3HU probe (Endocavity type probe)



Biopsy needle guides for Mindray C6-2GU probe (Microconvex type probe)



Biopsy needle guides for Mindray L12-3E probe (Linear type probe)



Biopsy needle guides for Mindray 6LB7s probe (Transperineal type probe)

Accurate, Safe and Efficient Operation



· Attach



· Angle Adjustment



· Needle Insertion



· Detach

List of Biopsy Needle Guides for Mindray

MODEL	Transducer	Transducer Type	Gauge Size/Guide Channel Length	MODEL	Transducer	Transducer Type	Gauge Size/Guide Channel Length
JSM-015	C5-2, C5-2s, C5-2E	Convex	11-23G/3.5cm	JSM-446	SL10-3U	Linear	11-23G/3.5cm
JSM-006	3C5A, 3C5, 3C5s, 3C5P	Convex	11-23G/3.5cm	JSM-447	LM18-4WU	Linear	11-23G/3.5cm
JSM-001	35C50EA, 35C50EB, 35C50P	Convex	11-23G/3.5cm	JSM-074	3C1, 3C1s	Microconvex	11-23G/3.5cm
JSM-045	SC5-1E	Convex	11-23G/3.5cm	JSM-075	65C15EA, 6C2, 6C2s, 6C2P	Microconvex	16G/18G/3.2cm
JSM-022	C5-1(s,E,U), SC6-1(s,E,U), C6-2(s,P), SC5-1N(s), SC5-1s	Convex	11-23G/3.5cm	JSM-123	C4-1, C4-1U, C4-1s	Microconvex	14-22G/3.5cm
JSM-219	C7-3E	Convex	11-23G/3.5cm	JSM-118	C11-3(s, E, U)	Microconvex	11-23G/3.5cm
JSM-031	SC5-1U, SC5-1E	Convex	11-23G/3.5cm	JSM-024	C6-2GE, C6-2Gs, C6-2GU	Microconvex	11-23G/3.5cm
JSM-265	SC8-2U, SC8-2s	Convex	11-23G/3.5cm	JSM-298	35C20EA	Microconvex	11-23G/3.5cm
JSM-348	C6-1	Convex	11-23G/3.5cm	JSM-448	SC6-1GU	Microconvex	11-23G/3.5cm
JSM-444	SCM7-1U	Convex	11-23G/3.5cm	JSM-229	C9-3	Microconvex	16G/18G/3.2cm
JSM-445	SC9-2U	Convex	11-23G/3.5cm	JSM-221	P7-3, P7-3s, P7-3E, P7-3P, P7-3U	Phased Array	11-23G/3.5cm
JSM-116	C6-2	Convex	16G/18G/3.2cm	JSM-443	SPM6-1U	Phased Array	11-23G/3.5cm
JSM-007	7L4(s,A,P), 7L4B(P), 10L4(s), 7L5(s,P), L7-3(s,E), L11-4(s), L12-4(s), L12-3E, L13-3(s), L14-6N(s,E), L14-6W(s,E,U), 75L53EA(HA)	Linear	11-23G/3.5cm	JSM-011	2P2(s), 2P2P, P4-2(s), P4-2E, P4-2NE, P4-2Ns, SP5-1s, SP5-1E, SP5-1U, SP5-1N, SP5-1Ns	Phased Array	11-23G/3.5cm
JSM-088	75L38EA, 75L38EB, 75L38P	Linear	11-23G/3.5cm	JSM-109	6LB7(s), 6LB7E, 6LB7P, 6LE7(s), 6LE7P, 65EL60EA	Transperineal	18G/3cm
JSM-036	LM15-6U	Linear	11-23G/3.5cm	JSM-246	ELC13-4U, ELC13-4s	Transperineal	18G/3cm
JSM-127	L14-5W, L14-5Ws, L14-5WE, L14-5WU	Linear	11-23G/3.5cm	JEM-028	65EC10(EA,EB,EC,ED,HA), 65EB10EA, 6CV1(s,P), CB10-4(s,P,E,A), V10-4(s), V10-4B(s,P), V11-3(E,s), V11-3B(E), V11-3W(s,E), E9-3, V11-3H(E,U,s)	Endocavity	16-18G/14.5cm
JSM-237	LM14-6E, LM14-6s, LM16-4U	Linear	11-23G/3.5cm	JEM-063	DE10-3(s,E,U), DE10-3W(s,E,U)	Endocavity	16-18G/19cm
JSM-284	L10-3E, L10-3s	Linear	11-23G/3.5cm	JEM-043	V11-3HB	Endocavity	16-18G/20.5cm
JSM-214	L14-3U, L13-3Ns	Linear	11-23G/3.5cm	JEM-064	V11-3H(E,U,s)	Endocavity	16-18G/21cm
JSM-126	L11-3U	Linear	11-23G/3.5cm	JEM-008	DE11-3U, DE11-3P, DE11-3s, DE11-3E, DE11-3	Endocavity	16-18G/18cm
JSM-287	L20-5, L20-5s, L20-5E, L20-5U	Linear	11-23G/3.5cm	JEM-092	SV10-2U	Endocavity	16-18G/15cm
JSM-204	L16-4Hs, L16-4HE, L16-4HU	Linear	16G/18G/3.2cm	JEM-080	DE11-3WU, DE11-3WE, DE11-3Ws	Endocavity	16-18G/18cm
JSM-264	L9-3U, L9-3E, L9-3s	Linear	11-23G/3.5cm	JEM-042	E9-4	Endocavity	16-18G/14.5cm
JSM-282	7LT4(s), 7LT4P, 75LT38EA, 75LT40HA	Linear	18G/20G/3.2cm	JSM-281	D8-2E, D8-2U	4D	14-22G/3.5cm
JSM-274	L14-6(s), L14-6P, L14-6E, 10L24EA	Linear	14-22G/3.5cm	JSM-144	D8-4U	4D	16G/18G/3.2cm
JSM-074	L11-3VNs	Linear	11-23G/3.5cm				
JSM-146	L12-3RC, L12-3RCs	Linear	20-23G/1cm				
JSM-305	L14-3WU, L14-3WE, L14-3Ws, L15-3WU, L13-3WE, L13-3WU	Linear	11-23G/3.5cm				

(Please contact InnoFine for more models.)