

Biopsy Needle Guides for SIUI



- Integral design without any accessories
- High grade medical stainless steel
- Artificial polishing
- Supporting accurate and angle adjustable biopsy
- Matching a wide range of needle gauges
- Simple pressing, quick release



About Us

As a China national high-tech enterprise, InnoFine are dedicated to R&D, manufacturing and sales ultrasound interventional products and consumables medical products since its establishment in 2013, providing safe and reliable, cost-effective products and overall solutions for thousands of hospitals in more than 100 countries around the world, products covering Ultrasound Needle Guides, Probe Cover, Gel, Biopsy Needle etc.

InnoFine will continue to adhere to the spirit of "Be innovative and forge ahead, Strive for excellence and perfection, Customer oriented ", and make unremitting efforts to promote the development of global interventional ultrasound field.

Cooperation Ultrasound Brand (Part)

| | | | | | |
|--|--|--|--|--|--|
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

Reference Sites (Part)

| | | | | |
|---|--|---|--|---|
| Mount Elizabeth Novena Hospital, Singapore | Apollo Hospitals, India | Semmelweis University Medical Center, Hungary | Kinderwunschzentrum Heidelberg, Germany | Orchid Fertility, UAE |
| Chinese PLA General Hospital | Peking Union Medical College Hospital | West China Hospital Sichuan University | Shengjing Hospital China Medical University | Tongji Hospital, Huazhong University of Science and Technology |
| Shanghai Ninth People's Hospital | Xiangya Hospital Central South University | The First Affiliated Hospital Sun Yat-sen University | Shenzhen People's Hospital | Shenzhen Hospital Peking University |

Classial Product (Part)



Biopsy needle guides for SIUI L8LC(50) probe (Linear type probe)



Biopsy needle guides for SIUI V6L11C probe (Endocavity type probe)



Biopsy needle guides for SIUI C3L60 probe (Convex type probe)



Biopsy needle guides for SIUI C3L60 probe (Convex type probe)



Biopsy needle guides for SIUI L8L38 probe (Linear type probe)



Biopsy needle guides for SIUI L8L50 probe (Linear type probe)

Accurate, Safe and Efficient Operation



· Attach



· Angle Adjustment



· Needle Insertion



· Detach

List of Biopsy Needle Guides for SIUI

| MODEL | Transducer | Transducer Type | Gauge Size/Guide Channel Length |
|---------|---|-----------------|---------------------------------|
| JSM-150 | C3L60, C3LC, C3L60C | Convex | 11-23G/3.5cm |
| JSM-336 | C3LC, C3L60C | Convex | 11-23G/3.5cm |
| JSM-336 | C3L60 | Convex | 11-23G/3.5cm |
| JSM-141 | L8L50 | Linear | 14-22G/3.5cm |
| JSM-319 | L8L38, L10LC | Linear | 14-22G/3.5cm |
| JSM-222 | L8LC(50), L8LN(50) | Linear | 11-23G/3.5cm |
| JSM-278 | L8L38C | Linear | 18-20G/3cm |
| JEM-002 | V6L11C, V6L11K, V6I11F, V5H11, V6LA, V6LC | Endocavity | 16-18G/15cm |
| JEM-090 | V6LC (4D) | Endocavity | 16-18G/18cm |

(Please contact InnoFine for more models.)