

Biopsy Needle Guides for Wisonic



- Integral design without any accessories
- High grade medical stainless steel
- Artificial polishing
- Supporting accurate and angle adjustable biopsy
- Matching a wide range of needle gauges
- Simple pressing, quick release



About Us

As a China national high-tech enterprise, InnoFine are dedicated to R&D, manufacturing and sales ultrasound interventional products and consumables medical products since its establishment in 2013, providing safe and reliable, cost-effective products and overall solutions for thousands of hospitals in more than 100 countries around the world, products covering Ultrasound Needle Guides, Probe Cover, Gel, Biopsy Needle etc.

InnoFine will continue to adhere to the spirit of "Be innovative and forge ahead, Strive for excellence and perfection, Customer oriented ", and make unremitting efforts to promote the development of global interventional ultrasound field.

Cooperation Ultrasound Brand (Part)

Reference Sites (Part)

 Mount Elizabeth Novena Hospital, Singapore	 Apollo Hospitals, India	 Semmelweis University Medical Center, Hungary	 Kinderwunschzentrum Heidelberg, Germany	 Orchid Fertility, UAE
 Chinese PLA General Hospital	 Peking Union Medical College Hospital	 West China Hospital Sichuan University	 Shengjing Hospital China Medical University	 Tongji Hospital, Huazhong University of Science and Technology
 Shanghai Ninth People's Hospital	 Xiangya Hospital Central South University	 The First Affiliated Hospital Sun Yat-sen University	 Shenzhen People's Hospital	 Shenzhen Hospital Peking University

Classial Product (Part)



Biopsy needle guides for Wisonic S4-2 probe (Phased Array type probe)



Biopsy needle guides for Wisonic C5-1B probe (Convex type probe)




Biopsy needle guides for Wisonic EV10-4 probe (Endocavity type probe)



Biopsy needle guides for Wisonic L15-4 probe (Linear type probe)



Biopsy needle guides for Wisonic L15-4B probe (Linear type probe)



Biopsy needle guides for Wisonic L10-5 probe (Linear type probe)

Accurate, Safe and Efficient Operation



· Attach



· Angle Adjustment



· Needle Insertion



· Detach

List of Biopsy Needle Guides for Wisonic

MODEL	Transducer	Transducer Type	Gauge Size/Guide Channel Length
JSM-295	C6-1B-H	Convex	11-23G/3.5cm
JSM-059	C7-3B, C5-1B, C5-2B	Convex	11-23G/3.5cm
JSM-050	C5-1, C5-2, C7-3	Convex	11-23G/3.5cm
JSM-427	CS6-2-U	Convex	11-23G/3.5cm
JSM-428	C9-2G-U	Convex	11-23G/3.5cm
JSM-296	L15-4WB-H	Linear	11-23G/3.5cm
JSM-297	L15-6-H	Linear	11-23G/3.5cm
JSM-040	L15-4B, L15-4NB	Linear	11-23G/3.5cm
JSM-041	L15-4	Linear	11-23G/3.5cm
JSM-032	L10-5	Linear	11-23G/3.5cm
JSM-426	L15-4W-U	Linear	11-23G/3.5cm
JSM-453	LH15-6-H	Linear	16-18G/1cm
JSM-299	C7-2	Microconvex	11-23G/3.5cm
JSM-420	C11-3-U	Microconvex	11-23G/3.5cm
JSM-452	HU-11C3	Microconvex	11-23G/3.5cm
JSM-483	C8-3	Microconvex	11-23G/3.5cm
JSM-026	S4-2	Phased Array	11-23G/3.5cm
JSM-020	P4-1	Phased Array	11-23G/3.5cm
JEM-028	EV10-4, EV11-3-H	Endocavity	16-18G/14.5cm
JEM-071	EV11-2-U	Endocavity	16-18G/15cm

(Please contact InnoFine for more models.)